DECLARACIÓN DE CONFORMIDAD CE
CE DECLARATION OF CONFORMITY
DÉCLARATION DE CONFORMITÉ CE

El fabricante | The manufacturer | Le fournisseur:

AUXILAB, S.L.

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Thanks for purchasing this equipment. We wish you make the most of its possibilities. We recommend you care this equipment following this manual.

Nahita develops its product according to directives of CE and emphasizing the ergonomics and safety.

Quality of materials used in manufacturing and a proper used will let you enjoy the equipment for many years.

An improper use of the equipment could cause accidents, electric shocks, short circuits, fires, injuries, etc.

TO GET THE BEST RESULTS AND A LONG LIFE PLEASE READ CAREFULLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT

Take into account the following:

- ◆ This manual is an inseparable part of nahita digital pH meter 903 so it should be available to all users.
- ◆ It should be handled with care, avoiding sudden movements, shocks or manipulation with sharp objects. Avoid spilling liquid on it.
- ◆ Never dismantle the equipment, do not try to repair it by yourself, you will lose the warranty and maybe you could suffer damages.
- ◆ If you have some doubts, please contact your distributor; installation, setting up or operation. You can also send your doubts or suggestions to NAHITA technical service: asistencia@auxilab.es or calling 807117040 (0.35 €/min.).
- ◆ This equipment is covered by law of warranties and consumer goods (10/2003).
- ◆ Fuses or accessories, as well as the loss of them, and wear pieces are not covered by the warranty.
- ◆ Manipulation by unauthorized personnel causes total loss of warranty.
- ◆ Make sure you keep the purchase invoice to have the right to claim or assistance. In case of sending the equipment to the technical service please enclose a copy of the invoice.
- ◆ The manufacturer reserves the right to make changes and improvements on this manual and on the equipment.

ATTENTION! WE DO NOT ADMIT EQUIPMENTS TO REPAIR IF THEY ARE NOT CLEAN AND DISINFECTED

10. ACCESSORIES

pH electrodes. Connector BNC

<i>Code</i>	<i>Diaphragm</i>
Glass body	
90900651	Ceramic
90900652	Ceramic
90900953	With ground joints
90900955	With ground joints
Plastic body	
90900218	Fiber
90900219	Fiber
90900900	PTFE Ring
90900901	PTFE Ring

pH buffers.

<i>Code</i>	<i>Description</i>
90900000	pH: 4.01, 7.00, 10.01

Appendix I. Diaphragms types of reference electrode / Applications

<i>Diaphragm</i>	<i>Flow KCl sat 25°C</i>	<i>Application</i>
Ceramic type Z	14 µl / h	Aqueous solutions with conductivity Z > 10 µS/ cm
Ceramic type D	2.5...5 µl / h	Aqueous solutions with conductivity Z > 10 µS/ cm
Glass with ground joints	0.5...2 ml / h	
Porous glass (dust of glass)	1.5...4ml / h	For organic environments
PTFE Ring	-	Dirty or oily environments

INSTRUCTIONS ABOUT ENVIRONMENTAL PROTECTION

Do not deposit the equipment in regular trash when its life cycle had finished; take it to a collection point for recycling electrical and electronic equipment.

Does not contain dangerous or toxic elements to humans but an inappropriate elimination will damage the environment.

Materials are recyclable. Recycling you are making a very important contribution to environment.

Please contact the administration of your community to know collection points.

ATTENTION! WE DO NOT ADMIT EQUIPMENTS TO REPAIR IF THEY ARE NOT WELL DISINFECTED.

9. TROUBLE SHOOTING

<i>PROBLEM</i>	<i>CAUSE</i>	<i>SOLUTION</i>
Digits are irregular	Wrong connection of the electrode	Check the connections
Irregularity in the measurements	Dirty or damaged electrode	Clean electrode Replace electrode
Display shows "Err" when on	Anomaly during check	Turn off the equipment, wait for a while and turn on again. If the problem persists, please, contact the Technical Service
Display shows "Err" in mV value and right temperature value	Input potential exceeds measurement range	Insert protector device in the connection of the electrode to re-establish the system
Button "pH/mV" does not let select measurement pH mode	When measurement pH range exceed ± 20.00 display indicates "Err" and continues to mV measurement mode	Reduce input signal Insert protector device to the electrode connection to re-establish the system

INDEX OF LANGUAGES

Spanish 2-16
 English 17-30
 French..... 31-45

INDEX OF CONTENTS

1.FEATURES OF THE EQUIPMENT..... 18
 2.THEORETICAL FUNDAMENTALS..... 19
 3.DESCRPTION 20
 4.TECHNICAL SPECIFICATIONS 21
 5.INSTALLATION / SETTING UP 21
 6.CALIBRATION 23
 7.Ph-METER FUNCTIONS 26
 8.MAINTENANCE AND CLEANING..... 27
 9.TROUBLE SHOOTING 29
 10.ACCESSORIES..... 30
 APPENDIX I: Diaphragm types of referent electrode / Features..... 30
 CE CERTIFICATE..... 46

1. FEATURES OF THE EQUIPMENT

Digital pH-meter is indicated to make measurings of pH at compensating temperature and electrochemical potential (mV) of watery solutions, so we can get diverse information, for example:

- ◆ Direct measurement of concentrations of anions and cations.
- ◆ Direct measurement of dissolved gases as NH₃, O₂, CO₂
- ◆ Acidity and alkalinity in water
- ◆ Calcium in water
- ◆ pH in wter, slurries, suspensions
- ◆ Proteins in water

Its features extend to many fields of work, we emphasized the following:

- ◆ Environment: land, sludge, sediments, waters.
- ◆ Biotechnology: fish farms, marine farmings, hydroponics
- ◆ Food: milks, sauces, beverages, meat, canned
- ◆ Pharmaceutical / cosmetic: gels and shampoos, deodorants, creams
- ◆ Petrochemical and general chemistry

2. THEORETICAL FUNDAMENTALS

The technique called direct potentiometry is the measure of the activity (or concentration) of a chemical species, measuring directly the potential with the one which is connected, through a known logarithmic function called Nernst equation.

$$E = E_0 + 2,3 (RT/nF) \log a_i$$

E = measured potential (mV)

E_0 = Potential of reference (mV)

R = gas constant

T = absolute temperature (K)

n = valency of the ion

F = Faraday constant

a_i = ionic activity

- ◆ In this equation “ E ” is the developed potential by the system, “ E_0 ” is the potential of reference, “ R ” is the gas constant, “ T ” is the absolute temperature, “ n ” is the number of electrons transferred and “ F ” is the Faraday constant.
- ◆ a_i represents the activity of the ion of interest. Ionic activity is related to the concentration (c_i) by the activity factor (f_i).

$$a_i = f_i c_i$$

- ◆ at constant temperature, Nernst equation is:

$$E = E_0 + S a_i$$

“ S ” is the slope of the electrode

From this equation we can deduce that making measurements of potential respect of a reference electrode, we can know the activity and, therefore, the ion concentration. The most known application of direct potentiometries is the use of what is known as Selected Ions Electrodes (SIE).

To obtain valid analytical potentiometry measurements, one of the electrodes has to be constant potential and it has not to suffer changes between experiments. Electrode that satisfies this condition is known as reference electrode. Due to the stability of the reference electrode, any change in the potential of the system will be thanks to the contribution of the other electrode, called working or indicator electrode.

Reference electrode can be an individual electrode or be incorporated into the indicator electrode (combined electrode). See Appendix I

- ◆ pH-meter has to be calibrated regularly in addition to the conditions indicated in “Calibration” section, by calibration buffers.
- ◆ “Life expectancy” of a pH electrode is about one year, this period varies depending on the quality of the electrode and its use: according to number of daily measurements, measurement at the extremes of pH scale and temperature, and especially, the maintenance.
- ◆ If duration is too short, may not be using the appropriate electrode!
- ◆ Additional liquid of combined reference electrode is a solution KCl 3 mol/L, that it can be added by the hole on the upper electrode if it is refillable.
- ◆ Electrodes must be kept into its protective storage with the corresponding electrolyte (avoiding accidental breakages of glass body electrodes), keeping this way the bulb (the most sensitive part of the electrode) always wet. We recommend KCl 4 M solution.
- ◆ Regeneration of the electrode when it detects a decrease of the slope after a prolonged use:
 - Submerge the bult in HF 4% for 3~5s
 - Clean it with distilled water
 - Submerge it in HCl 0.1 M
- ◆ Do not leave the electrode submerged in distilled water, protein or hydrofluoric acid solution (HF) for so long.
- ◆ The diaphragm that carries the electrode has not to be exposed to water-insoluble contaminant materials (organic oils, silicone, resins, pigments, etc), they can blocked the diaphragm and prevent proper contact among solutions, being the main reason errors in measurements.
- ◆ When pH-meter is not being used, always place the protective device output BNC to prevent that, in any case, the circuit is opened and damage the equipment.
- ◆ Inspect the electrode and the wire. Wire used to connect to the pH-meter must be intact, without cracks, like the electrode. If not, replace the electrode. Always use original parts and components. Others may be similar but their use can damage the equipment.
- ◆ Never dismantle the equipment, do not try to repair it by yourself, you will lose the warranty and maybe you could suffer damages.
- ◆ In case of failure, contact your supplier for repairing through Nahita Technical Service.

Cleaning

Important: each time you make a measurement you should clean the electrode with distilled water, avoiding this way that samples suffer contamination of previous solution.

- ◆ Electrode has to be perfectly clean with distilled water. Do not use a cloth (because it can be charge electrostatically), use a paper free of fuzz (filter paper, cellulose).
- ◆ For cleaning different pieces of the equipment, we recommend use neutral products that they do not contain abrasives agents, such as a few drops of soap for cleaning glass lowered in water.
- ◆ Do not spill the prepared on the equipment, you should wet a soft cloth with the product.

Possible causes of error:

- ◆ Alkaline error: ordinary glass electrodes become sensitive to alkaline materials with pH value higher than 9.
- ◆ Acid error: with a real pH value less than zero, obtained values used to be higher.
- ◆ Dehydration: false results. You can not make the measurements.
- ◆ Temperature: pH measurement changes with the temperature, it can be corrected though system of internal temperature compensation equipment.

Safety

- ◆ pH-meters have to be used by qualified people, they should know the equipment and its use.
- ◆ Do not place the equipment in areas near heat sources (burners, torches...) or expose directly to sunlight, etc.
- ◆ When pH-meter is not being used, put always protector device in the output BNC (1.3) to avoid that circuit is opened and damages to electric components. This way you also avoid dust and humidity.
- ◆ BNC output should be kept clean and dry.
- ◆ In order to protect the pH-meter and guarantee correct measurements, always check the electrode/s connections on the back before making the measurement.
- ◆ Once finished working we recommend to unplug the equipment because electric components in the adapter has a small power consumption which means premature aging of equipment, possible accidents due to heat production as well as an unnecessary waste of energy, etc.
- ◆ For any cleaning manipulation, checking or replacement of components (ex. Fuse replacement) is essential shutdown and unplugged the equipment.
- ◆ Do not try to repair it by yourself, you will lose the warranty and cause damages in the operation as well as injuries to people (burns, injuries...) and damages to the wiring.
- ◆ Do not let water in the control panel, even if it is adequately insulated. If for any reason you suspect that there is water or any liquids in it, unplugged the equipment.
- ◆ pH-meter has a Schuko wire; it have to be connected to grounded current.
- ◆ Manufactured according to European Directives about electric safety, electromagnetic compatibility and safety in machines.


8. MAINTENANCE AND CLEANING

Every previously mentioned use rules have no value without an appropriate maintenance:


- ◆ Follow instructions and warnings given in this manual.
- ◆ Keep this manual to be consulted.

3. DESCRIPTION

- | | |
|---|------------------------------------|
| 1.1 LCD Screen | 1.4 Reference electrode connection |
| 1.2 Keyboard | 1.5 Fuse |
| 1.3 BNC connection of indicator electrode / Protective device | 1.6 Switch on / off |
| | 1.7 Current |

**Front panel**

- | | |
|---|--|
| 2.1 pH / mV: selecting measurement mode (pH/potential) | 2.5 ENTER: Confirmation of parameters / system reset |
| 2.2 TEMP: temperature selection | |
| 2.3 SLOPE: adjustment pH (2 ^o calibration point) | |
| 2.4 STD: adjustment pH (1 ^o calibration point) | |

4. TECHNICAL SPECIFICATIONS

Code	50903000
pH range	0~14.00 pH
mV range	0~±1999 mV
Resolution	0.01 pH / 1 mV
Accuracy of the electronic unit	±0.01 pH/ ±1 mV ±1 bit
Accuracy of the measurement (Detector)	±0.02 pH ± bit
Repetitive accuracy (Electronic unit)	pH: 0.01/ mV: 1 mV
Repetitive accuracy (Detector)	≤ 0.01 pH
Temperature compensation range	0-60 °C
Temperature compensation accuracy	±0.01 pH
Input current	≤ 1x10 ⁻¹² A
Input impedance	≥ 1x10 ¹² Ω
Stability of the electronics unit	±0.01 pH ± 1 bit /3 h
Supply	AC (220±22) V / (50±1) Hz
Fuse	0.5 V
Dimensions	290x210x95 mm
Weight	1.5 Kg

Among its main features we can emphasize:

- ◆ Wide LCD screen, easy to read and understand
- ◆ Measurement functions of pH at compensated temperature and electrical potential (mV) with automatic polarity indicator
- ◆ Solution temperature selector
- ◆ Automatic detection of standard buffer solution (4.00, 6.86 y 9.18 pH)

5. INSTALLATION / SETTING UP

Preliminary inspection

- ◆ Unpacked the pH-meter, remove the plastic and the polyspan protection.
- ◆ Be sure the pH-meter does not present transport damages, if so, please contact immediately with your distributor in order to make the claim within the prescribed period by carrier.

Returns of equipments may be made before 15 days after the sending and come back complete, in their original packing with all accessories and documents included

7. pH-METER FUNCTIONS

Important: each time you make a measurement you should wash the electrode with distilled water, avoiding samples suffer contamination of previous solution.

Attention; do not leave the electrode submerged in distilled water

pH / mV measurement

Electrodes known as pH electrodes are systems for electrochemical potential measuring (mV) of hydrogen ion (H⁺) allowing us to know their activity or concentration (pH) in determined aqueous solution.

Proceed as follows:

- ◆ Remove protector device of connection BNC (1.3) and connect the combined electrode of pH.
- ◆ Measure the temperature of the problem solution with a thermometer by button "TEMP" (2.2) adjust the screen temperature to the temperature of the problem solution.
- ◆ Select measurement mode by button "pH/mV" (2.1).
- ◆ Submerge the electrode in the problem solution, shake it slowly and proceed to the reading of the pH measurements and potential on the display.

Note: reference electrode may be an individual electrode or be incorporated to indicator electrode (combined electrode). In case of being individual, both must be plugged separately, first to the reference electrode connection (1.4) and second to BNC connection (1.3).

Potential electrode measurement of selective ion.

When it is required to determine the concentration of a specific ion different from H⁺, it is necessary to use Selective Electrodes of Ions (SEI), their membrane responds more or less to a particular ion.

Process as follows:

- ◆ Remove protective device from BNC connection (1.3) and connect selective electrode of ion.
- ◆ Remove electrode case and clean it with distilled water and after with problem solution.
- ◆ Select potential measurement mode mV by button "pH/mV" (2.1).
- ◆ Submerge the electrode in problem solution, shake it softly and make the reading of the measurement. Polarity is showed automatically.

Note: if the measurement is out of detected range by pH-meter, display will indicate "Err".

(2.2) adjust the temperature on the screen and temperature of the buffer solution.

- ◆ Press button "SLOPE" (2.3) (indicative "Std yes" on the screen) and after press "ENTER" (2.5). pH-meter automatically detects present standard solution and shows its pH value and selected temperature. For example:

If you use a different buffer solution, adjust by button SLOPE (▲/▼) (2.3), value pH on the screen to pH solution value and working temperature.

Note: To know this pH, look up Temperature-pH table that it is included in the back of the envelopes of the buffer.

- ◆ Press button "ENTER" to finish the calibration. Screen indicates "Std ok" and shows the value of the slope and the potential E0 that are automatically registered. Finally it returns to measurement pH mode.
- ◆ Wash the electrode with distilled water and dry it with filter paper or cellulose.

Note: after calibration, the system records new reference data and automatically deletes the previous calibration data, the equipment is ready for use.

pH-meters have to be calibrate when:

- pH-meter is new
- Electrode is new
- Long periods of time with no use
- Long use periods, in this case calibrate it every 24 hours
- After using chemical aggressive substances (acids, alkalis, etc.)

Reset slope of the electrode

When errors occur in the slope of the system owing to various reasons (electric power cut, mistakes during calibration, etc) it is possible to re-establish its value:

- ◆ Re-calibrating the equipment following previous paragraph
- ◆ By button "ENTER"
 - In measurement mode (pH/mV) press "ENTER" (2.5) for approx. 5 seconds, display shows the indication "SYS rst" (fisxd / flashing).
 - Press again "ENTER" (2.5) and displays shows default values of the slope and E0 corresponding to reference points pH=4.00 and pH=9.18 (Slope: 100.0 %, E0: 0 mV).
 - Finally it returns to potential measurement mode ("0mV").

- ◆ Check the accessories you have to receive with the equipment:

- Buffer solution pH 4.00, 6.86 and 9.18
- Standar Schuko wire
- User's manual
- Warranty certificate

Installation

Before you start using your pH-meter, you should familiarize yourself with its components and basic fundamentals, as well as, the functions of its controls.

READ THIS MANUAL BEFORE WORKING WITH THE EQUIPMENT IN ORDER TO OBTAIN THE MAXIMUM FEATURES AND A HIGHER LIFE

- ◆ Place the pH-meter on an horizontal and stable table, leaving space enough around the working area.
- ◆ Do not place the equipment in areas near heat sources (burners, torches...) or expose directly to sunlight, etc.
- ◆ Keep the pH-meter outside interference from magnetic fields.
- ◆ Insert the Schuko power cord alternating current (AC) to the current base 220 V, 50 Hz +/- 10% supplied with grounded, and at the other end to the power connector of the pH-meter.

Neither the manufacturer nor the distributor assume any responsibility for damages to the equipment, installation or injuries to people owing to the inobservance of the correct connection procedure. Tension must be 220 V, 50 Hz ± 10%.

- ◆ When the pH-meter is not being used, always put the protective device in the output BNC (1.3) to avoid, the circuit is opened and damage the equipment.
- ◆ When you are not going to use the pH-meter for long periods of time, make sure it is disconnected from the power supply.

Setting up

- ◆ Connect the wire to the supply current (1.7).
- ◆ Remove the protector device of BNC connection (1.3) and connect the combined electrode.
- ◆ Press the power switch (1.6).

Screen shows the following indications:

Being:

“S=100.0”; value of the slope in percent.

“E0=0.0 mV”: value of electric potential at the isopotential point system.


Both values are defined after the calibration of equipment.

Finally it comes into potential measurement mode and shows “0 mV” and the last programmed temperature.

- ◆ The equipment is ready for use.

Note: we recommend waiting 30 minutes. for stabilization of the pH-meter before proceeding to make measurements.

System operates following this diagram:

6. CALIBRATION

Calibration can be made by:

- ◆ Automatic detection of standard buffer solutions (4.00 pH, 6.86 pH y 9.18 pH) or manual adjustment with different buffer solutions.
- ◆ One or two points of calibration.

Note: clean very well the electrode before its calibration because resulting measurements, based on the calibration of a dirty electrode will be wrong. Be sure the used buffer sample for the standardization of the pH-meter is in good conditions.

Calibration at a point

It is used when it is not necessary to make precise measurements.

- ◆ Prepare the standard buffer solution pH=6.86
- ◆ Select measurement mode of pH by button pH/mV (2.1)
- ◆ Clean the electrode with distilled water and submerge in the solution
- ◆ Measure the temperature of the buffer solution with a thermometer using button “TEMP” (2.2) adjust the temperature on the screen with the temperature of the buffer solution (ej. 25°C).
- ◆ Press button “STD” (2.4) (indicative “Std yes” on the screen) and after “ENTER” (2.5). The pH-meter detects automatically the present standard solution and it shows its pH value and selected temperature.

If you used a different buffer solution, please adjust by button STD (▲/▼) (2.4), the pH value on the screen to the value of pH solution to the working temperature.

Note: to know this pH, look up Temperature-pH table that it is included in the back of the envelopes of the buffer.

- ◆ Press button “ENTER” (2.5); screen indicates “Std ok” and shows the slope value and potential E0 are automatically recorded until next calibration. Finally it returns to pH measurement mode.
- ◆ Wash the electrode with distilled water and dry it with filter paper or cellulose.

Note: after calibration, the system records new reference data and automatically deletes the previous calibration data, the equipment is tested for its use.

The slope for the calibration at a point is considered by default 100.0%

Calibration at two points

Generally it is used to calibrate the slope of the electrode

- ◆ Prepare standard buffer solutions pH=4.00 and pH=9.18.
- ◆ Proceed to make the calibration with the solution pH=4.00 following steps 2. to 5. of “Calibration at one point”.
- ◆ Press button “ENTER” (2.5); pH-meter returns to pH measurement mode.
- ◆ Wash the electrode with distilled water and dry it with filter paper or cellulose.
- ◆ Submerge the electrode in the solution pH=9.18.
- ◆ Measure the temperature of buffer solution with a thermometer using “TEMP” button