Crash Cart Equipment



Updated March 2006: D. Tucker, RPh, BCPS with K. Wahby, PharmD and R. Simpson, PharmD

Objectives

- To familiarize the pharmacist with some of the non-medication contents of the crash cart drawers
- To describe the uses of some of the drawer components

Outside of the Cart

- Monitor and paddles
 - Defibrillation
 - Evaluate heart rate
- Suction equipment
- Needle box
- Examination gloves
- Oxygen tank

Monitor and Paddles









Suction Apparatus





HWH

HUH/RIM

Needle Boxes



HWH



HUH/RIM

Oxygen Canister





Airway Equipment

- Intubation equipment
- Endotracheal tubes
 - If cuffed can use with ventilator
- CO₂ detector
 - To make sure tube is in airway


Intubation Equipment

Curved blade

Extra batteries

Handle


Straight blade



Endotracheal Tubes & Airway



CO₂ Detector



Breathing Equipment

- Face mask
- Manual resuscitator
- Oxygen equipment for monitoring and delivery

Face Mask



Adult Manual Resuscitator



Oxygen Flowmeter




Nasal Oxygen Cannula



Circulation Equipment

- To insert a central line
- To start peripheral lines
- To test blood
 - Venous
 - ABGs (put on ice and deliver to lab STAT)

Central Venous Catheter (CVC) Insertion Kit



IV Bags and Tubing



Tegaderm Transparent Dressing



ABG Kit



Butterfly



Blood Collection Vials



Syringes and Alcohol Swabs




Safety Needles



Personal Equipment

- Face mask
- Cover gown
- Gloves
 - Sterile
 - Examination

Face Mask



Sterile Gloves



Miscellaneous Equipment

- Scissors
- Hemostat
- Irrigation equipment
- Sponges
- Extension cord

Bandage Scissors



Hemostat



NS and Irrigation Syringe



Sponges



Extension Cord

